

# life - POINT



CV. MITRA PERKASA  
Jl. Kapi Janula 15 D / 14, Malang  
Telp: (0341) 5473579, 7039809  
Fax: (0341) 722839

## BIPHASIC DEFIBRILLATOR

CE  
1984

# PRO BIPHASIC 270 JOULE

## Technical Specifications

Wave Form	Biphasic
Operating Mode	Synchronous / Asynchronous
Energy Selection	2 - 270 Joule
ECG Derivations	DI-DII-DIII-AvR-AvL-AvF-C1...C6
ECG Filters	Mon1-Mon2-Strong-Surgery-Diagnostic
Heart Rate Alarm Limits	15 - 300
Freeze	Yes
Monitor Type	LCD
Monitor Resolution	320x240
Monitor Size	5.5 inch
Monitor Speed	25mm/sec
Printer Type	Thermal Printer
Printer Speed	25-50mm/sec
Paper Size	50mm x 30m
Max. Discharging	100+
ECG Monitoring Time	4 Hours
Working Conditions	0/40 °C ve %30-95Rh
Storing Conditions	-20 / +55 °C
Dimensions	403x152x324 mm
Weight	8,5 Kg
Class	II B
Main Supply Voltage	AC230V/1A 50-60 Hz
Direct Working From Main Voltage	Optional
Pulseoximeter	Optional
Standart Accessories	ECG Cable and Paper ECG Electrode - Power Cable

## Opsiyonel Aksesuarlar



Stand



Trolley